FR

Pansement PolyMem® Surgical™

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EN

PolyMem® Surgical™ Dressing

Step 1. Remove one side of adhesive release liner.

Step 2. Apply over wound. Slowly remove second half of release liner while pressing adhesive to the skin.

Step 3a. Pinch dressing slightly where two white strips come together at edge of dressing. This will cause cover sheet to start to lift.

Step 3b. Remove cover sheets one at a time. Gently smooth thin film as each cover sheet is being removed.

DESCRIPTION:

PolyMem polymeric membrane dressings are made of moisture-loving polyurethane matrix that contains (1) safe wound cleanser, (2) glycerin moisturizer, and (3) superabsorbents. Dressings have a thin polyurethane backing which breathes and helps keep out contamination. These dressings may be used as a combined primary and secondary dressing or as a secondary dressing. This highly conformable and absorbent dressing in MAX thickness draws wound fluid, which is known to contain natural growth factors and nutrients, to the wound site.

PolyMem Silver products incorporate very fine silver particles and will not stain skin. Silver acts inside the dressing to protect from microbial contamination. Approximate silver content by area 207 µg/cm²

Drug-free PolyMem dressings help to :
• effectively manage and heal wounds
• help to reduce risk of surgical site infection
• absorb fluid and provide a moist healing environment
• minimize and relieve persistent and procedural wound pain
• help to accelerate post-surgical recovery and rehabilitation
• reduce edema, bruising, and the spread of inflammation into surrounding undamaged tissues

INDICATIONS:

Under the direction of a healthcare professional, for the management of partial and full thickness surgical wounds healing by primary intention.

PRECAUTIONS:

- For external use only • PolyMem dressings are not compatible with oxidizing agents (hydrogen peroxide or hypochlorite solutions)
- Do not re-sterilize • Do not use if package is open or damaged
- Single use • Do not apply the dressing under tension • Choose a dressing that ensures coverage of the incision by the PolyMem pad; the adhesive should not come in contact with the incision
- PolyMem dressings are not intended to replace skin closures such as staples, sutures or tissue adhesive • Allow tissue adhesive to dry thoroughly before applying, or pad may become attached
- Prior to performing imaging procedures, determine if the dressing should be removed • Avoid contact with electrodes or conductive gels • Topical treatments are not recommended in conjunction with PolyMem dressings • Be alert for signs of infection or maceration • Do not use and discontinue use on people who show signs of sensitivity, irritation, or allergy from the dressing or its materials
- Intended for one use on one person. Do not reuse. Reuse may degrade properties, spread contamination, and increase risk of infection.

PolyMem dressings may be used when signs of infection are present if proper medical treatment that addresses the cause of the infection has been implemented.

Store at room temperature ≤30°C (86°F). Excursions permitted ≤5°C (41°F) but should be minimized. **Dispose** as typical for biohazardous waste. **Dimensions** are approximate.

EN

DIRECTIONS:

Initial Application:

1. Prepare the incision site or chronic wound according to protocol or as directed by a physician or other ordering clinician. Dry periwound completely to ensure long-lasting dressing adhesion.
2. Choose a dressing that ensures direct contact of the PolyMem Pad with the entire incision or chronic wound so that the membrane pad is ¼" – 2" (0.6 – 5 cm) larger than the incision or open wound.
3. Best results are achieved if the PolyMem pad covers all inflamed or damaged areas surrounding the incision or the chronic wound. Do not apply under tension.
4. For dry chronic wounds, moisten the wound or the dressing slightly.
5. Outline the chronic wound margins on the outside of the dressing to help determine when to change.
6. For knee incision, knee should be at approximately a 45 to 60 degree angle for optimal application to reduce risk that leg flexion / extension will result in skin damage.

Use:

1. For chronic wounds, a dramatic increase in wound fluid may be observed during the first few days due to the fluid attraction. This is not uncommon and indicates that the dressing is working.
2. Keep the dressing dry and in place when bathing. Change the dressing if it gets wet.

Dressing change indications:

1. **For incisions,** perform initial dressing change 3-5 days after primary closure. Other reasons to change:
 - The pad is 75% saturated by fluid or blood
 - The patient or wound has signs of infection
 - At discharge, follow guidelines of facility or surgeon
 - If the incision is not completely epithelialized, replace with another PolyMem dressing. The wound usually does not require additional cleansing before application.

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- Continued application after epithelialization improves long-term outcomes. Apply new dressings weekly until all signs of inflammation have resolved. Replace if the dressings become wet.
- 2. For an exuding chronic wound, fluid will become visible through the top. Change before fluid reaches the approximate wound margin, or when good practice dictates, or after no more than 7 days. Change immediately if immune response is observed.
- 3. More frequent changes may be needed due to a compromised immune system, diabetes, infection, hypergranulation, macerated tissue, or when desiring to facilitate the removal of non-viable tissue.

Dressing Change:

1. Dressing Change process: Lift the edge of the dressing. Using the 180-degree removal technique peel the dressing back onto itself. Pressing the skin at the skin/adhesive junction helps assure gentle removal. Continue throughout the length and width of the dressing.
2. If coagulated blood or dried exudate has caused the dressing to adhere to skin or hair, use a small amount of saline or water to slowly gently soften and loosen. Clean the intact skin around the wound.
3. Manage dressing changes according to protocol or as directed by a physician or other ordering clinician. In most cases, when using PolyMem, there is no need to disturb or cleanse an incision site during changes unless the incision site is infected or contaminated. PolyMem is designed to continuously clean the incision site and does not leave residue that needs to be removed.

Excessive cleaning may injure regenerating tissue and delay healing.

- During normal healing of chronic wounds, the wound may appear to become larger during the first few dressing changes because the body is naturally removing non-viable tissue. If the wound site continues to become larger, consult an expert.
- PolyMem dressings help support new blood vessel formation so it is not uncommon to see blood-stained wound fluid and dressings during changes.
- Be alert for signs and symptoms, such as, increased redness, inflammation, pain, warmth, odor, whitening, fragility, excessive tissue formation, or lack of healing. Consult a professional for problem situations, such as, infection, irritation, sensitivity, allergy, maceration, or hypergranulation.

Any serious incident that has occurred in relation to this device should be reported, by the user and/or patient, to Ferris Mfg. Corp. and the competent authority of the Member State in which the user and/or patient is established.